

BIOLOC™ CELL CULTURE BAG CONFIGURATIONS



| VOLUME CAPACITY | SERIES DESIGN | CCA FOR ADHERENT CELLS | CCN FOR SUSPENSION CELLS | SINGLE PORT | DUAL PORT | 3+ PORT |
|-----------------|---------------|------------------------------|--------------------------------|----------------|--------------|------------|
| 3mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | | |
| | 3 | | | | | |
| | 4 | ✓ | ✓ | ✓ | | |
| 10mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | | |
| | 3 | | | | | |
| | 4 | ✓ | ✓ | ✓ | | |
| 30mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | |
| | 3 | | | | | |
| | 4 | ✓ | ✓ | ✓ | ✓ | |
| 32mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | |
| | 3 | | | | | |
| | 4 | ✓ | ✓ | ✓ | ✓ | |
| 70mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | |
| | 3 | ✓ | ✓ | ✓ | ✓ | |
| | 4 | ✓ | ✓ | ✓ | ✓ | |
| 72mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | |
| | 3 | ✓ | ✓ | ✓ | ✓ | |
| | 4 | ✓ | ✓ | ✓ | ✓ | |
| 120mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | |
| | 3 | ✓ | ✓ | ✓ | ✓ | |
| | 4 | ✓ | ✓ | ✓ | ✓ | |
| 198mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | |
| | 3 | ✓ | ✓ | ✓ | ✓ | |
| | 4 | ✓ | ✓ | ✓ | ✓ | |
| 240mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 3 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 4 | ✓ | ✓ | ✓ | ✓ | ✓ |
| 250mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 3 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 4 | ✓ | ✓ | ✓ | ✓ | ✓ |
| 275mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 3 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 4 | ✓ | ✓ | ✓ | ✓ | ✓ |
| 290mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 3 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 4 | ✓ | ✓ | ✓ | ✓ | ✓ |
| 300mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 3 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 4 | ✓ | ✓ | ✓ | ✓ | ✓ |
| 325mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 3 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 4 | ✓ | ✓ | ✓ | ✓ | ✓ |
| 500mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 3 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 4 | ✓ | ✓ | ✓ | ✓ | ✓ |
| 750mL | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 3 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 4 | ✓ | ✓ | ✓ | ✓ | ✓ |
| 1L | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 3 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 4 | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2L | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 3 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 4 | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2.5L | 1 | ✓ | ✓ | | ✓ | |
| | 2 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 3 | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 4 | ✓ | ✓ | ✓ | ✓ | ✓ |